KORONA[®]


DELEILA

- D** **DIAGNOSE-WAAGE** Gebrauchsanleitung
- GB** **DIAGNOSTIC SCALE** Instruction for Use
- FR** **PÈSE-PERSONNE IMPÉDANCEMÈTRE** Mode d'emploi
- ES** **BÁSCULA PARA DIAGNÓSTICO** Instrucciones para el uso
- IT** **BILANCIA DIAGNOSTICA** Istruzioni per l'uso
- GR** **ΔΙΑΓΝΩΣΤΙΚΗ ΖΥΓΑΡΙΑ** Οδηγίες χρήσεως
- RUS** **ΔΙΑΓΝΩΣΤΙΚΕΣ ΒΕΣΕΣ** Инструкция по применению
- PL** **WAGA DIAGNOSTYCZNA** Instrukcja obsługi
- NL** **PERSONENWEEGSCHAAL MET**
EXTRA FUNCTIES Gebruiksaanwijzing

 Hans Dinslage GmbH
Riedlinger Straße 28
88524 Uttenweiler, Germany
Tel.: +49 7374 / 915766
www.korona-net.de

[D] Deutsch

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Dear Customer

You have chosen a quality product from Korona. Please read these instructions for use carefully and keep them for later use, be sure to make them accessible to other users and observe the information they contain. This diagnostic scale is designed to make a significant contribution to your health.

With kind regards
Your Korona team

1. Interesting facts

The measuring principle of the diagnostic scale

For weight measurement and bone mass measurement, the results are shown in 100 g increments (0.2 lb, 1/4 st.). The measuring results of the body fat, body water and muscle percentage are shown in 0.1% increments. With this process, an imperceptible, completely safe and harmless current enables the determination of body tissue proportions within seconds. The body fat percentage and other physical relationships in the body can be determined by measuring the electrical resistance (impedance), and the inclusion in a calculation of constants and individual parameters (age, height, gender, degree of activity). Muscle tissue and water have good electrical conductivity, and therefore a lower resistance. Bones and fatty tissue, on the other hand, have low conductivity, as fat cells and bones conduct hardly any current as a result of their very high resistance. Please be aware that values obtained from the diagnostic scale represent only an approximation of actual analytical medical data. Only a specialist physician can accurately determine body fat, body water, muscle percentage and bone structure using medical procedures (e.g. computed tomography).

Weighing correctly

If possible, always weigh yourself at the same time of day (ideally in the morning), after going to the toilet, on an empty stomach and without clothes to achieve comparable results.

Only the long-term trend counts.

Short-term weight fluctuations within a few days are usually only due to a loss of liquid. However, body water plays an important role in your well-being. The interpretation of the results depends on the changes in the:

1. total weight
2. percentage of fat
3. percentage of water
4. percentage of muscle

and on the length of time during which these changes take place. Rapid changes over several days should be distinguished from mid-term changes (over several weeks) and long-term changes (months). As a basic rule, SHORT-TERM weight changes almost exclusively represent changes in water content, while mid-term and long-term changes can ALSO concern the percentage of fat and muscle.

- If you experience a short-term drop in weight, however the percentage of body fat increases or remains the same, you have only lost water, e.g. after a training session, going to the sauna or following a diet limited to rapid weight loss.
- If you experience a mid-term weight increase, while the percentage of body fat drops or remains the same, you might have built up your muscle mass or your body may have retained more water (e.g. as the result of taking medication or a weaker cardiac function in elderly people).
- If your weight and percentage of body fat drop simultaneously, your diet is working and you are losing fat mass. Ideally, you can support your diet with physical exercise, fitness or power training. This can result in a mid-term increase in the percentage of muscle in your body.

And remember: Physical exercise is the basis for a healthy body.

2. Important notes – keep for later use!

2.1 Safety notes

- **The scale must not be used by persons with medical implants (e.g. heart pacemakers). Otherwise their function could be impaired.**
- Do not use during pregnancy.
- Caution! Do not stand on the scale with wet feet and do not step on the scale when its surface is wet – danger of slipping!
- Keep children away from packing materials (danger of suffocation).
- Batteries are highly dangerous if swallowed. Keep batteries and scale out of reach of small children. If batteries are swallowed, get medical help immediately.
- Batteries must not be charged or reactivated by other means, taken apart, thrown into a fire or short-circuited.

2.2 General notes


- The unit is for personal use only and is not intended for medical or commercial applications.
- Please note that technically related measuring tolerances are possible, as the scale is not calibrated for professional medical applications.
- The diagnostic scale can be preset for an age range of 10–100 years and heights of 100–220 cm (3' 03"–7' 02"). The scale can accept a maximum load of 150 kg (330 lb, 23 st.). The weight measurement is displayed in increments of 100 g (0.2 lb, 1/4 st.) angezeigt. For weight measurement and bone mass measurement, the results are shown in 100 g increments (0.2 lb, 1/4 st.). The measuring results of the body fat, body water and muscle percentage are shown in 0.1% increments.
- On delivery, the scale is set to "cm" and "kg". At the back of the scales is a switch with which you can switch between "inches", "pounds" and "stones" (lb, st.).
- Place the scale on a firm level floor; a firm floor covering is required for correct measurement.
- Repairs may only be carried out by customer service or authorized dealers. Before submitting any complaint, first check the batteries and replace them if necessary.
- All scales comply with EC Directive 2004/108/EC + supplements. Should you have any questions concerning the use of our scale, please contact your dealer or customer service.
- Please dispose of the blanket in accordance with the directive 2002/96/EG – WEEE (Waste Electrical and Electronic Equipment). If you have any queries, please refer to the local authorities responsible for waste disposal.

3. Operation

3.1 Batteries

If present, pull the battery insulating strip off the battery compartment cover or remove the battery's protective film and insert the battery according to the polarity. If the scale fails to operate, remove the battery completely and insert it again. Your scale is equipped with a "replace battery" indicator. When the scale is operated with a battery which is too weak, "LU" appears in the display field and the scale automatically switches off. The battery must be replaced in this case (4 x 1,5 V, AAA, Mikro).

The used, completely drained standard and rechargeable batteries must be disposed of in specially marked collection containers, at toxic waste collection points or electrical product retailers. You are legally obligated to dispose of the batteries.

Note: Batteries containing pollutants are marked with the following symbols:

Pb = Battery contains lead, Cd = Battery contains cadmium, Hg = Battery contains mercury.

Place the scale on a firm, flat surface; a firm floor covering is a basic requirement for correct weight measurement.

3.2 Measuring weight only

Audibly, briefly and firmly tap the scale with your foot to activate the vibration sensor! The display “8888” appears as a self test (Fig. 1), until “0.0” (Fig. 2) is displayed. Now the scale is ready to measure your weight. Stand on the scale without moving about and distribute your weight evenly on both legs. The scale immediately begins to measure your weight. Your measured weight is displayed shortly after this (Fig. 3). The scale switches off a few seconds after you step down from it.

Fig. 1

Fig. 2

Fig. 3

Fig. 4

Fig. 5

Fig. 6

Fig. 7

Fig. 8

3.3 Measurement of body weight, as well as body fat, body water and muscle percentage, bone mass and Kcal

First your personal parameters must be stored;

- Tap your foot quickly, forcefully and audibly on the scale to activate the vibration sensor and wait until the display “0.0” (Fig. 2) appears. Then press “SET”. Now the first memory preset flashes in the display (Fig. 4). Select the desired preset by pressing the buttons “up” ▲ or “down” ▼ and confirm the desired preset with the “SET” button.
- The sex setting “Woman” or “Man” appears on the display (Fig. 5). Select your sex by pressing the buttons “up” ▲ or “down” ▼ and confirm by pressing “SET”.
- “Age” in years then flashes in the display (Fig. 6). Set the desired value by briefly pressing or holding the buttons “up” ▲ or “down” ▼ and confirm it by pressing “SET”.
- Now height in cm (inch) flashes in the display (Fig. 7). Set the desired value by briefly pressing or holding the buttons “up” ▲ or “down” ▼ and confirm it by pressing “SET”.
- Now the setting for your individual level of activity or fitness appears on the display (Fig. 8). Select the desired level of activity by pressing the buttons “up” ▲ or “down” ▼ and confirm the desired level of activity with the “SET” button.

Explanation of the 5 levels of activity

- **Activity level 1:** No physical activity while overweight
- **Activity level 2:** Little physical activity while overweight
- **Activity level 3:** No or little physical activity
(Less than 20 min of physical exercise twice a week, e.g. going for a walk, light gardening, gymnastic exercises)
This mode should also be selected before the start of a training or dietary program. After approx. 6–10 weeks you can switch to the next level of activity.
- **Activity level 4:** Moderate physical activity
(20 minutes of physical exercise 2 to 5 times a week, e.g. jogging, cycling, tennis etc.)
This mode should be selected as soon as you notice a general improvement in your well-being/fitness in the course of a training/dietary program. After 8–12 weeks of continuous training/diet, you can change to the next level of activity.
- **Activity level 5:** Intensive physical activity
(Daily intensive training or physical work, e.g. intensive running, construction work etc.) For this mode you should have an overall high level of mobility, endurance and power and maintain this level over a longer period of time.

The scale has 10 memory presets for basic data with which the users, e.g. each member of the family, can store and call up their personal settings at any time.

After all parameters have been entered, you can now measure your weight, body fat and the other data.

- **Audibly, briefly and firmly tap the scale with your foot to activate the vibration sensor!**
- Then press “SET”. Now the first memory preset flashes in the display (Fig. 4).
- Select the memory preset where your basic personal data are stored by pressing the ▲ or ▼ button repeatedly. The data are displayed consecutively until “0.0” (Fig. 2) appears.
- If the memory space does not flash, press “SET” again.
- Step onto the scale barefoot and make sure you’re standing on both electrodes. First your body weight is determined and displayed.

- Remain on the scale and avoid moving; your body fat and water, and your muscle and bone mass are now being analysed (Fig. 9). This may take a few seconds.

Fig. 9

Important: Your feet, legs, calves and thighs must not touch each other, as otherwise the measurement cannot be performed correctly.

The following data are displayed:

- **FAT** Body fat percentage, in %
- **TBW** Water percentage, in %
-  Muscle percentage, in %
-  Bone mass, in kg
- **KCAL** Kcal

Now all measured values are displayed consecutively and the scale switches off.

3.4 Tips for using the scale

Important for calculating body fat, body water and muscle percentage, as well as bone mass:

- The measurement may only be carried out while barefoot and it is helpful if the soles of your feet are slightly damp. Completely dry soles can result in unsatisfactory results, as they have insufficient conductivity.
- Stand still during the measurement.
- Wait several hours (6–8) after unusually strenuous activity.
- Wait approx. 15 minutes after getting out of bed so that the water in your body can be distributed.

The measurement is not reliable for:

- Children under approx. 10 years of age.
- Persons with fever, undergoing dialysis, with symptoms of edema or osteoporosis.
- Persons taking cardiovascular medication. Persons taking vasodilating or vasoconstricting medications.
- Persons with substantial anatomical deviations in the legs relative to their total height (leg length considerably shorter or longer than usual).

Body fat guide

The following body fat levels provide you with a guideline (for further information, please consult your doctor!). A lower level is often found in athletes. Depending on the type of sports, the intensity of training and the person's physical constitution, levels can be achieved that are even lower than the specified guidelines.

Women

Age	Underfat	Healthy	Overfat	Obese
10–39	<21%	21–25%	26–31%	>32%
40–55	<22%	22–26%	27–32%	>33%
56–100	<23%	23–27%	28–33%	>34%

Men

Age	Underfat	Healthy	Overfat	Obese
10–39	<18%	18–23%	24–29%	>30%
40–55	<19%	19–24%	25–30%	>31%
56–100	<20%	20–25%	26–31%	>32%

According to the World Health Organization (WHO 2001), the percentage of body water should lie within the following ranges!

Women: 50–55%

Men: 60–65%

Children: 65–75%

In persons with a high percentage of body fat, the percentage of body water very frequently lies below the specified guidelines.

As the percentage of muscle differs from individual to individual, there are no generally valid guidelines.

Bone mass

Like the rest of our body, our bones are subject to the natural development, degeneration and ageing processes. Bone mass increases rapidly in childhood and reaches its maximum between 30 and 40 years of age. Bone mass reduces slightly with increasing age. You can reduce this degeneration somewhat with healthy nutrition (particularly calcium and vitamin D) and regular exercise.

With appropriate muscle building, you can also strengthen your bone structure.

Note that this scale will not show you the calcium content of your bones, but will measure the weight of all bone constituents (organic substances, inorganic substances and water).

Attention: Please do not confuse bone mass with bone density. Bone density can be determined only by means of a medical examination (e.g. computer tomography, ultrasound). It is therefore not possible to draw conclusions concerning changes to the bones and bone hardness (e.g. osteoporosis) using this scale. Little influence can be exerted on bone mass, but it will vary slightly within the influencing factors (weight, height, age, gender).

AMR

The active metabolic rate (AMR) is the amount of energy required daily by the body in its active state. The energy consumption of a human being rises with increasing physical activity and is measured on the diagnostic scale in relation to the degree of activity entered (1–5).

To maintain your existing weight, the amount of energy used must be reintroduced into the body in the form of food and drink. If less energy is introduced than is used over a longer period of time, your body will obtain the difference largely from the amount of fat stored and your weight will reduce.

If, on the other hand, over a longer period of time more energy is introduced than the total active metabolic rate (AMR) calculated, your body will be unable to burn off the excess energy, and the excess will be stored in the body as fat and your weight will increase.

4. Incorrect measurement

If the scale detects an error during weighing, “Err” or “Err2” appears in the display.

If you step onto the scale before “0.0” appears in the display, the scale will not operate properly.

Possible causes of errors:	Remedy:
– The maximum load-bearing capacity of 150 kg (330 lbs) was exceeded.	– Only weigh the maximum permissible weight.
– The electrical resistance between the electrodes and the soles of your feet is too high (e.g. with heavily callused skin).	– Repeat weighing barefoot. Slightly moisten the soles of your feet if necessary. Remove the calluses on the soles of your feet if necessary.
– Your body fat lies outside the measurable range (less than 5% or greater than 50%).	– Repeat weighing barefoot. – Slightly moisten the soles of your feet if necessary.
– Your percentage of water lies outside the measurable range (less than 30% or greater than 90%)	– Repeat weighing barefoot. – Slightly moisten the soles of your feet if necessary.

5. Cleaning/care

The scale should be cleaned occasionally with a damp cloth. Do not use abrasive detergents and never immerse the scale in water.

Protect the scale against hard knocks, moisture, dust, chemicals, major temperature fluctuations and heat sources which are too close (stove, heating radiators).

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