KORONA[®]

Britta

D	GLASWAAGE Gebrauchsanleitung	2
GB	GLASS SCALE Instruction for Use	5
FR	PÉSE-PERSONNE EN VERRE Mode d'emploi	8
ES	BÁSCULA DE VIDRIO PARA PERSONAS Instrucciones para el uso	11
IT	BILANCIA IN VETRO Istruzioni per l'uso	14
GR	ΓΥΑΛΙΝΗ - ΥΓΑΡΙΑ Οδηγίες χρήσεως	17
RUS	СТЕКЛЯННЫЕ ВЕСЫ Инструкция по применению	20
PL	WAGA SZKLANA Instrukcja obsługi	23
NL	WEEGSCHAAL MET GLAZEN WEEGPLATEAU Gebruiksaanleiding	26
P	BALANÇA DE VIDRO Instruções de utilização	29
TR	CAM TERAZI Kullanma Talimatı	32
DK	GLASVÆGT Brugsanvisning	35
S	GLASVÄG Bruksanvisning	38
FIN	LASIVAAKA Käyttöohje	41
CZ	SKLEŇNÁ VÁHA Návod k použití	44
SLO	STEKLENA TEHTNICA Navodila za uporabo	47
H	ÜVEGMÉRLEG Használati útmutató	50

Korona electric GmbH
Hauptstraße 169, 59846 Sundern, Germany
Tel.: +49 2933 / 90284-80
www.korona-electric.de

Service:

Korona Service, Am Steinbach 9, 59872 Meschede-Enste

1 Anwendungsbedingungen

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1 Conditions of use

You should weigh yourself with as little clothes as possible. Make sure to weigh yourself at the same time of day each time. Since the body is subject to natural fluctuations (e.g. dehydration from sports or sauna, intake of food and fluids), constant conditions should be ensured when measuring, so as to receive comparable results.

KORONA electric GmbH shall not accept liability for damages or losses caused by the bathroom scales or for third party claims. This product is intended for private use only. It is not designed for professional use in hospitals or medical facilities.

2 First use

Please make sure to check the packaging and device for any visible damage before using the scales for the first time. The scales must not be used if any damage is detected. Please contact us immediately. Our customer service representatives will be happy to help.

Safety guidelines: This device is not to be used by persons (including children) with limited physical, sensory or mental abilities or who are lacking experience and understanding, unless they are supervised by persons, who are responsible for their safety or have been instructed in the use of the scales. It should be ensured that children do not play with the device.

Place the scales on a flat and solid surface (tiles, parquet, etc.). Incorrect measurements may occur on carpets.

Stand in the centre of the scale surface to prevent tipping,

Attention! There is a danger of slipping if the scales' surface is wet or if feet are wet.

Ensure that no objects fall onto the scales; falling objects could break the glass and damage the scales.

Use only a damp cloth to clean the scales, Never use solvents or abrasive cleaning agents, Do not immerse in water,

Preparation: Open the lids of the battery compartments on the underside of the scales and on the reverse of the portable display module to remove the battery safety strips, Ensure the correct polarity of the batteries (+/-).

Pressing the small black button UNIT inside the battery compartment of the display module will toggle the unit of weight kg <-> lb (kg = kilogram / lb = pound). Please ensure that the scales are configured to show the unit of weight that you normally use.

If you would like to mount the display module on the wall, remove the battery compartment and attach it to the wall via the relevant screw opening. You can then suspend the display module from this bracket.

3 Weighing function

IMPORTANT! The scales are operated via an IR receiver. Please remember that the portable display and the scales must be directly connected while using the scales. There must not be any obstacles between the scales and the display. The IR range between the scales and portable display is limited.

The scales are equipped with a Step-On function. Hold the display module in your hand and step

onto the scales. Stand still and do not hold on to anything. Ensure that there is no obstacle between the scales (IR emitter) and the display module (IR receiver). The display will initially show 88:8.8, followed by your actual weight after a few seconds. Once the scales have stabilised, your weight will be displayed 3 times in succession. Step off the scales. The scales will switch off automatically after a few seconds.

4 Error messages

Overload - The 180 kg maximum load capacity has been exceeded. Remove the load from the scales immediately to avoid damage to the load cells.

Batteries in the display module are low - Please replace the battery with an appropriate battery type (1 x CR2032).

Batteries in the scales' are low - Please replace the battery with an appropriate battery type (1 x CR2032).

5 Technical data

Load capacity x gradation: Max. 180 kg x 100 g
Dimensions: 302 mm x 302 mm x 17 mm
Batteries: 2 x 1 pce. CR2032 (supplied) Subject to technical charges

6 Warranty

KORONA electric GmbH agrees to rectify any material or manufacturing defects by offering free repair or replacement, for a period of 5 years from the purchase date. In case of a warranty claim, please return the scales to your stockist, including proof of purchase and stating the reason for the claim.

Note: Extreme sources of electromagnetic radiation, such as those caused by the operation of radio equipment in the vicinity of the scales, can affect the value displayed. The product will return to proper function as soon as the interference has been eliminated. Note: the unit may need to be restarted.
Legal information obligation regarding the disposal of batteries

Do not dispose of batteries in domestic waste. As a consumer, you are legally required to return used batteries. You may return your used batteries to public collection points in your community or to any retailer that sells batteries of the same type.

Note

Pb = Battery contains lead,
Cd = Battery contains cadmium
Hg = Battery contains mercury

CR (Li); AA (Al, Mn); AAA (Al, Mn)
Safety information for lithium batteries (CR2032/CR2430)

1. Lithium batteries are lithium-manganese cells. If you are not intending to use the scales for a long period of time, you should remove the batteries from the device.
2. Never open the lithium batteries or throw them into fire or expose them to shocks, as harmful vapours could escape or batteries might explode.
3. Never touch leaking batteries with your bare hands:

4. If your eyes or hands are contaminated, rinse with a lot of water; in case of skin or eye irritation, consult a doctor.
5. Only dispose of cells that are completely discharged or packaged protected against polarity reversal in accordance with local disposal regulations.
6. Do not subject the cells to direct sunlight and extreme heat; otherwise, there is a risk of overheating.
7. Maintain a storage temperature of $< 40^{\circ}\text{C}$.

Disposal of electrical and electronic devices

Follow local regulations regarding the disposal of materials. The appliance should be disposed of according to Regulation 2002/96/EC-WEEE (Waste Electrical and Electronic Equipment). In case of queries, contact the municipal authorities responsible for waste disposal in your area.

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